# ***Biotest MMB Listeria monocytogenes (R)***

## **Manual**

Version 06.09

**Biotest MMB Listeria  
monocytogenes (R)  
Art. No. 944 032**

für den real-time PCR-Nachweis von *L. monocytogenes*-DNA/  
for the real-time PCR detection of *L. monocytogenes* DNA

## Inhaltsverzeichnis / Contents

### DEUTSCH

<b>EINLEITUNG .....</b>	<b>3</b>
<b>NORMEN .....</b>	<b>3</b>
<b>KIT-VERSIONEN (V, R, LC) .....</b>	<b>3</b>
<b>SENSITIVITÄT .....</b>	<b>4</b>
<b>KOMPONENTEN UND LAGERUNG .....</b>	<b>4</b>
<b>QUALITÄTSKONTROLLE .....</b>	<b>5</b>
<b>HINWEISE UND VORSICHTSMAßNAHMEN .....</b>	<b>5</b>
<b>PROTOKOLL .....</b>	<b>5</b>
<b>ERGEBNISSE .....</b>	<b>7</b>
<b>WEITERE INFORMATIONEN .....</b>	<b>8</b>
<b>BESTELLINFORMATIONEN .....</b>	<b>9</b>
<b>ENGLISH</b>	
<b>INTRODUCTION .....</b>	<b>10</b>
<b>NORMS AND REGULATIONS .....</b>	<b>10</b>
<b>SELECTION OF KIT VERSION (V, R, LC) .....</b>	<b>10</b>
<b>SENSITIVITY .....</b>	<b>11</b>
<b>COMPONENTS AND STORAGE .....</b>	<b>11</b>
<b>QUALITY CONTROL .....</b>	<b>12</b>
<b>WARNING AND PRECAUTIONS .....</b>	<b>12</b>
<b>PROTOCOL .....</b>	<b>12</b>
<b>INTERPRETATION OF RESULTS .....</b>	<b>14</b>
<b>FURTHER INFORMATION .....</b>	<b>15</b>
<b>ORDER INFORMATION .....</b>	<b>16</b>

#### Trademarks and Licensing

Applied Biosystems, VIC are registered trademarks; StepOne, FAM, ROX are trademarks of Applied Biosystems Corporation or its subsidiaries in the US and/or certain other countries.

Bio-Rad, iCycler iQ registered trademark, iQ5 is a trademark of Bio-Rad, Inc.

BHQ is a registered trademark of Biosearch Technologies, Inc.

RotorGene is a registered trademark of Corbett research.

Eppendorf, Mastercycler are registered trademarks of Eppendorf AG.

Roche, LightCycler is a registered trademark of Roche Molecular Systems, Inc.

Stratagene, Mx3000P, Mx3005P, Mx4000P are registered trademarks of Stratagene.

All other trademarks are the property of their respective owners.

The application in which this product is used is covered by patents issued and applicable in the United States and certain other countries. Because purchase of this product does not include a license to **perform any patented application, users of this product may be required to obtain a patent license.**

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hycon@biotest.de

## Introduction

Listeriosis is a infectious disease with low incidence (in Germany 100-200 cases/year) but high mortality rates, which is caused by food contaminated with bacteria of the species *Listeria monocytogenes* (different serovars). Immunocompromised people, elderly, pregnant, foetuses and neonates are especially susceptible for Listeriosis.

The real-time PCR method enables, in contrast to classic microbiological culture methods, the rapid and sensitive testing of raw material and processed food for absence of *Listeria monocytogenes*.

The test procedure can be performed using all common real-time cyclers. Three distinct versions of the test kit are available for different instrument types (see chapter 3): The internal controls for exclusion of a potential inhibition of the PCR reactions are labelled with different suitable fluorescent dyes.

## Norms and Regulations

The Biotest MMB *Listeria monocytogenes* test kit fulfills the following norms and regulations:

ENORM EN ISO 11290-1 und ENORM EN ISO 11290-2 Ausgabe: 2005-02-01 - Microbiology of food and animal feeding stuffs - Horizontal method for the detection and enumeration of *Listeria monocytogenes* - Part 1 and part 2: Detection method (ISO 11290-1:1996 + Amd.1:2004)

DIN EN ISO 20837, 2006-08

Polymerase chain reaction (PCR) for the detection of pathogen microorganisms in food – Requirements to the preparation of samples for the qualitative detection (ISO 20837:2006) German version EN ISO 20837, 2006-08

DIN EN ISO 20838, 2006-08

Polymerase chain reaction (PCR) for the detection of pathogen microorganisms in food – Requirements to the amplification and detection of qualitative method (ISO 20838:2006); German version EN ISO 20838:2006)

DIN EN ISO 22174, 2005-05 = BVL L 00.00-45

Polymerase chain reaction (PCR) for the detection of pathogen microorganisms in food – General Requirements and Definitions (ISO 22174:2005); German version EN ISO 22174:2005

BVL L 00.00-95 (V), Technical rule, 2006-12 = §64 official research methods LFGB (= §35 LMBG)

Detection of food - qualitative proof of *Listeria monocytogenes* in food by - Polymerase chain reaction (PCR)

## Selection of Kit Version (V, R, LC)

All Biotest MMB *Listeria monocytogenes* kits contain FAM-labeled probes for the detection of *Listeria monocytogenes* DNA at 522 nm. The Biotest MMB *Listeria monocytogenes* (R) kit uses the fluorescent dye ROX for the detection of the internal amplification control. For use with other detection platforms different variants are available, containing fluorescent labels which are detected on other channels:

Biotest MMB <i>Listeria monocytogenes</i> (V):	VIC channel (553 nm)
Biotest MMB <i>Listeria monocytogenes</i> (R):	ROX channel (602 nm)
Biotest MMB <i>Listeria monocytogenes</i> (LC):	LightCycler detection (640 nm)

The following table provides information, which kit variant is recommended for the respective instruments:

Instrument Description	Kit version		
	(V)	(R)	(LC)
Applied Biosystems ABI 7000, 7300, 7500 7700, 7900 SDS, Step One	x		
BioRad iCycler iQ iQ5		x x	
Roche LightCycler 1.5 LightCycler 2.0 LightCycler 480		x	x x
Eppendorf Mastercycler ep realplex <sup>2</sup> Mastercycler ep realplex <sup>4</sup>	x x		
Stratagene Mx Series	x		
Corbett Research RotorGene / QIAGEN Rotor-Gene Q	x		

## Sensitivity

The detection limit is  $\leq 5$  DNA copies. The sensitivity of the whole procedure is influenced by the sample matrix, processing grade, DNA isolation and DNA content.

## Components and Storage

- 2x Listeria monocytogenes (R) Reaction Mix (1.0 ml)
- 1x Positive Control DNA (180  $\mu$ l)
- 1x FDE (Fluorescence Detection Enhancer) (140  $\mu$ l)
- 1x Taq DNA Polymerase (10  $\mu$ l)

Store all reagents at  $-20$  °C protected from light.

Required equipment and materials, not included in the kit:

- Real-time cycler with detection channels 522 nm and 602 nm
- Adjustable pipettes
- Sterile DNA-free filter tips
- Suitable real-time PCR consumables (plates, tubes, optical clear sealing options, capillaries)
- Centrifuge with suitable rotors for tubes, plates or capillaries (if necessary, with adapters).

## Quality Control

According to the Quality Management System, each lot of Biotest MMB *Listeria monocytogenes* Kits is tested against predetermined specifications to ensure consistent product quality.

## Warning and Precautions

For research use only, not for *in vitro* diagnostic use.

**Caution:** The test must be performed by well-trained and authorized laboratory technicians.

**Caution:** All reagents should be handled in accordance with good laboratory practice taking appropriate precautions. Do not pipette by mouth.

**Caution:** Do not use reagents past the expiration date printed on the label.

**Caution:** Do not use reagents with any evidence of turbidity or microbial contamination.

**Caution:** Post-PCR areas should be physically separated from Pre-PCR areas of the laboratory.

**Caution:** Pipettes used for Post-PCR manipulations should not be used for Pre-PCR manipulations.

## Protocol

### Isolation of *Listeria monocytogenes* DNA

For isolating *Listeria monocytogenes* DNA the Biotest MMB Prep *Listeria* Kit (Art. No. 944 515) is recommended.

## Amplification of the Target Sequence

### PCR Setup

For preparation of the master mix determine the number of samples and controls (positive control, negative control and extraction control). The master mix contains an internal amplification control to monitor inhibition of the PCR due to possible contaminants in the DNA sample.

It is recommended to prepare 10 % additional volume to the master mix in order to compensate pipetting losses.

Before use, all reagents must be thawed, vortexed and centrifuged briefly.

Example for calculation and preparation of 10 reactions:

Components for master mix	Volume per reaction	10 reactions (+10 %)
Listeria monocytogenes (R) Reaction Mix	18.8 µl	206.8 µl
FDE (Fluorescence Detection Enhancer)	1.1 µl	12.1 µl
Taq DNA Polymerase	0.1 µl	1.1 µl
<b>Total volume</b>	<b>20 µl</b>	<b>220 µl</b>

Mix thoroughly and spin down briefly.

Pipette 20 µl master mix to each reaction tube, well, or capillary.

Close the reaction vessel for the negative control directly.

Add 5 µl DNA sample. Close reaction vessel directly after sample transfer.

Add 5 µl Positive Control DNA to the last reaction vessel. Close reaction vessel directly after transfer.

Spin samples briefly to collect them at the bottom of the tube.

Place reaction vessels in the PCR instrument and start the cycling program.

## Cycler Set-up

For the Biotest MMB *Listeria monocytogenes* (R) kit use the following cycler set-up:

Initial Denaturation (HOLD)	5 min, 95 °C
Cycles	45
Denaturation	15 s, 95 °C
Annealing/Extension (CYCLE)	30 s, 60 °C
Temperature Transition Rate/ Ramp Rate	Maximum
Fluorescence Detection Setup	Detection: End of extension phase Reporter Dye 1 (522 nm): FAM (detection system) Reporter Dye 2 (602 nm): ROX (internal amplification control) Quencher Dye 1+2: BHQ  <b>Roche LC480</b> Reporter dye 1 (483, 533 nm) Reporter dye 2 (558, 610 nm)

## Interpretation of Results

The evaluation must be made according to the analysis program provided by the real-time PCR instrument manufacturer. All control reactions need to show the correct results.

A sample is evaluated **positive**, if the sample DNA shows amplification in the detection system (FAM).

A sample is evaluated **negative**, if the sample DNA shows no amplification in the detection system and the internal amplification control (inhibition control, ROX channel) of the sample is positive.

If the sample DNA and the internal amplification control are negative the sample contains a PCR inhibiting contaminant. Under these circumstances an evaluation of the samples is not possible. The DNA preparation for the sample must be repeated or improved.

## **Additional Information**

For any inquiries or to place an order, contact your local Biotest distributor.  
Domestic customers in Germany may contact Customer Service at:

Phone: +49 (0) 6103 801-496

Fax: +49 (0) 6103 801-505

Email: [mail@biotest.de](mailto:mail@biotest.de)

For technical or applicational inquiries please contact our Technical Support at:

Phone: +49 (0) 6221 72651-30

Fax: +49 (0) 6221 72651-39

Email [hycon@biotest.de](mailto:hycon@biotest.de)

For further questions please send an e-mail to [hycon@biotest.de](mailto:hycon@biotest.de).

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Biotest (Schweiz) AG
Biotest Italia S.r.l.
Biotest Diagnostics Corporation