# ***Biotest MMB Salmonella (V)***

## **Manual**

Version 06.09

**Biotest MMB Salmonella (V)**

**Art. No. 944 011**

für den real-time PCR-Nachweis von Salmonellen-DNA /  
for the real-time PCR detection of Salmonella DNA

## Inhaltsverzeichnis / Contents

### DEUTSCH

<b>EINLEITUNG .....</b>	<b>3</b>
<b>NORMEN .....</b>	<b>3</b>
<b>KIT-VERSIONEN (V, R, LC) .....</b>	<b>3</b>
<b>SENSITIVITÄT .....</b>	<b>4</b>
<b>KOMPONENTEN UND LAGERUNG .....</b>	<b>4</b>
<b>QUALITÄTSKONTROLLE .....</b>	<b>5</b>
<b>HINWEISE UND VORSICHTSMAßNAHMEN .....</b>	<b>5</b>
<b>PROTOKOLL .....</b>	<b>5</b>
<b>ERGEBNISSE .....</b>	<b>7</b>
<b>WEITERE INFORMATIONEN .....</b>	<b>7</b>
<b>BESTELLINFORMATIONEN .....</b>	<b>10</b>
<b>ENGLISH</b>	
<b>INTRODUCTION .....</b>	<b>10</b>
<b>NORMS AND REGULATIONS .....</b>	<b>10</b>
<b>SELECTION OF KIT VERSION (V, R, LC) .....</b>	<b>10</b>
<b>SENSITIVITY .....</b>	<b>11</b>
<b>COMPONENTS AND STORAGE .....</b>	<b>11</b>
<b>QUALITY CONTROL .....</b>	<b>12</b>
<b>WARNING AND PRECAUTIONS .....</b>	<b>12</b>
<b>PROTOCOL .....</b>	<b>12</b>
<b>INTERPRETATION OF RESULTS .....</b>	<b>14</b>
<b>FURTHER INFORMATION .....</b>	<b>15</b>
<b>ORDER INFORMATION .....</b>	<b>16</b>

#### Trademarks and Licensing

Applied Biosystems, VIC are registered trademarks; StepOne, FAM, ROX, TAMRA are trademarks of Applied Biosystems Corporation or its subsidiaries in the US and/or certain other countries.

Bio-Rad, iCycler iQ registered trademark, iQ5 is a trademark of Bio-Rad, Inc.

Rotorgene is a registered trademark of Corbett research.

Eppendorf, Mastercycler are registered trademarks of Eppendorf AG.

Roche, LightCycler is a registered trademark of Roche Molecular Systems, Inc.

Stratagene, Mx3000P, Mx3005P, Mx4000P are registered trademarks of Stratagene.

All other trademarks are the property of their respective owners.

The application in which this product is used is covered by patents issued and applicable in the United States and certain other countries. Because purchase of this product does not include a license to **perform any patented application, users of this product may be required to obtain a patent license.**

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[hycon@biotest.de](mailto:hycon@biotest.de)

## Introduction

Salmonellosis is a food-borne infection with one billion incidences annually worldwide (2005). This infectious gastric disease is caused by several strains of the bacterial species *Salmonella enterica* (serovar). The consumption of contaminated food derived from animals or plants are the main cause for Salmonellosis.

The real-time PCR method enables, in contrast to classic microbiological culture methods, the rapid and sensitive testing of raw material and processed food for absence of *Salmonella*.

The test procedure can be performed using all common real-time cyclers. Three distinct versions of the test kit are available for different instrument types (see chapter 3): The internal controls for exclusion of a potential inhibition of the PCR reactions are labelled with different suitable fluorescent dyes.

## Norms and Regulations

The Biotest MMB Salmonella test kit fulfills the following norms and regulations:

DIN EN ISO 20837, 2006-08

Microbiology of food and animal feeding stuffs - Polymerase chain reaction (PCR) for the detection of food-borne pathogens - Requirements for sample preparation for qualitative detection

DIN EN ISO 20838, 2006-08

Microbiology of food and animal feeding stuffs – polymerase chain reaction (PCR) for the detection of food-borne pathogens – Requirements for amplification and detection for qualitative methods

BVL L 00.00-98, Technical rule, 2007-04 (§64 LFGB)

Detection of food – qualitative proof of Salmonella in food – with polymerase chain reaction (PCR)

## Selection of Kit Version (V, R, LC)

All Biotest MMB Salmonella kits contain FAM-labeled probes for the detection of Salmonella DNA at 522 nm. The Biotest MMB Salmonella (V) kit uses the fluorescent dye VIC for the detection of the internal amplification control. For use with other detection platforms different variants are available, containing fluorescent labels which are detected on other channels:

Biotest MMB Salmonella (V):	VIC channel (553 nm)
Biotest MMB Salmonella (R):	ROX channel (602 nm)
Biotest MMB Salmonella (LC):	LightCycler detection (640 nm)

The following table provides information, which kit variant is recommended for the respective instruments:

Instrument Description	Kit-version		
	(V)	(R)	(LC)
Applied Biosystems ABI 7000, 7300, 7500 7700, 7900 SDS, Step One	x		
BioRad iCycler iQ iQ5		x x	
Roche LightCycler 1.5 LightCycler 2.0 LightCycler 480		x	x x
Eppendorf Mastercycler ep realplex <sup>2</sup> Mastercycler ep realplex <sup>4</sup>	x x		
Stratagene Mx Series	x		
Corbett Research RotorGene / QIAGEN Rotor-Gene Q	x		

## Sensitivity

The detection limit is  $\leq 5$  DNA copies. The sensitivity of the whole procedure is influenced by the sample matrix, processing grade, DNA isolation and DNA content.

## Components and Storage

- 2x Salmonella (V) Reaction Mix (1.0 ml)
- 1x Positive Control DNA (180  $\mu$ l)
- 1x Taq DNA Polymerase (10  $\mu$ l)

Store all reagents at  $-20$  °C protected from light.

Required equipment and materials, not included in the kit:

- Real-time cycler with detection channels 522 nm and 553 nm
- Adjustable pipettes
- Sterile DNA-free filter tips
- Suitable real-time PCR consumables (plates, tubes, optical clear sealing options, capillaries)
- Centrifuge with suitable rotors for tubes, plates or capillaries (if necessary, with adapters).

## Quality Control

According to the Quality Management System, each lot of Biotest MMB Salmonella Kits is tested against predetermined specifications to ensure consistent product quality.

## Warning and Precautions

For research use only, not for *in vitro* diagnostic use.

**Caution:** The test must be performed by well-trained and authorized laboratory technicians.

**Caution:** All reagents should be handled in accordance with good laboratory practice taking appropriate precautions. Do not pipette by mouth.

**Caution:** Do not use reagents past the expiration date printed on the label.

**Caution:** Do not use reagents with any evidence of turbidity or microbial contamination.

**Caution:** Post-PCR areas should be physically separated from Pre-PCR areas of the laboratory.

**Caution:** Pipettes used for Post-PCR manipulations should not be used for Pre-PCR manipulations.

## Protocol

### Isolation of Salmonella DNA

For isolating Salmonella DNA the Biotest MMB Prep Salmonella Kit (Art. No. 944 505) is recommended.

In sterile samples or samples with a very low microbial activity (e.g. heat pasteurized egg white powder) a contamination with not viable or not culturable *Salmonella* cells might be observed. In case of other conditions the samples show a fast metabolization of these cells.

To validate results in matrices with low or zero microbial activity, it is recommended to compare the Ct-value for the signal for the samples from at the start and at the end of the pre-enrichment (dynamic).

## Amplification of the Target Sequence

### PCR Setup

For preparation of the master mix determine the number of samples and controls (positive control, negative control and extraction control). The master mix contains an internal amplification control to monitor inhibition of the PCR due to possible contaminants in the DNA sample.

It is recommended to prepare 10 % additional volume to the master mix in order to compensate pipetting losses.

Before use, all reagents must be thawed, vortexed and centrifuged briefly.

Example for calculation and preparation of 10 reactions:

Components for master mix	Volume per reaction	10 reactions (+10 %)
Salmonella (V) Reaction Mix	19.9 µl	218.9 µl
Taq DNA Polymerase	0.1 µl	1.1 µl
<b>Total volume</b>	<b>20 µl</b>	<b>220 µl</b>

Mix thoroughly and spin down briefly.

Pipette 20 µl master mix to each reaction tube, well, or capillary.

Close the reaction vessel for the negative control directly.

Add 5 µl DNA sample. Close reaction vessel directly after sample transfer.

Add 5 µl Positive Control DNA to the last reaction vessel. Close reaction vessel directly after transfer.

Spin samples briefly to collect them at the bottom of the tube.

Place reaction vessels in the PCR instrument and start the cycling program.

## Cycler Set-up

For the Biotest MMB Salmonella (V) kit use the following cycler set-up:

Initial Denaturation (HOLD)	5 min, 95 °C
Cycles	45
Denaturation	15 sec, 95 °C
Annealing/Extension (CYCLE)	30 sec, 60 °C
Temperature Transition Rate/ Ramp Rate	Maximum
Fluorescence Detection Setup	Detection: End of extension phase Reporter Dye 1 (522 nm): FAM (detection gene) Reporter Dye 2 (553 nm): VIC (internal amplification control) Quencher Dye 1+2: TAMRA

## Interpretation of Results

The evaluation must be made according to the analysis program provided by the real-time PCR instrument manufacturer. All control reactions need to show the correct results.

A sample is evaluated **positive**, if the sample DNA shows amplification in the detection system (FAM).

A sample is evaluated **negative**, if the sample DNA shows no amplification in the detection system and the internal amplification control (inhibition control, VIC channel) of the sample is positive.

If the sample DNA and the internal amplification control are negative the sample contains a PCR inhibiting contaminant. Under these circumstances an evaluation of the samples is not possible. The DNA preparation for the sample must be repeated or improved.

## **Additional Information**

For any inquiries or to place an order, contact your local Biotest distributor.  
Domestic customers in Germany may contact Customer Service at:

Phone: +49 (0) 6103 801-496

Fax: +49 (0) 6103 801-505

Email: [mail@biotest.de](mailto:mail@biotest.de)

For technical or applicational inquiries please contact our Technical Support at:

Phone: +49 (0) 6221 72651-30

Fax: +49 (0) 6221 72651-39

Email [hycon@biotest.de](mailto:hycon@biotest.de)

For further questions please send an e-mail to [hycon@biotest.de](mailto:hycon@biotest.de).

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Biotest (Schweiz) AG
Biotest Italia S.r.l.
Biotest Diagnostics Corporation